


# HIGH PERFORMANCE LIGHTWEIGHT 37g

THE NEW GEUDER®  
ERGO-LIGHT ULTRASONIC  
HANDPIECE



MADE FOR MEGATRON® S3 SYSTEM

NEW ERGONOMIC DESIGN    INNER IRRIGATION    EXTREMELY LIGHT - ONLY 37 g



ORIGINAL SIZE

**G-30440**  
**ERGO-LIGHT**  
**ULTRASONIC HANDPIECE**  
for MEGATRON® S3 system

diameter 13 mm, length 125 mm, weight 37 g

## ADVANTAGES

- extremely light and thin titanium housing
- new ergonomic design with excellent grip
- consistent stroke over the whole lifetime
- inner irrigation tube
- 6 piezolayers for perfect efficiency
- for all phaco techniques including CMP and MICS

**Geuder**®  
Precision made in Germany

GEUDER AG  
Hertzstrasse 4  
69126 Heidelberg  
Germany

Phone: +49 6221 3066  
Fax: +49 6221 303122  
info@geuder.de  
www.geuder.de