

THE NEW HIGH-SPEED VITREOUS CUTTER MACH2

THE DOUBLE-BLADE EFFICIENCY.

**IMMEDIATELY
AVAILABLE**

**FAST CORE VITRECTOMY AND
SAFE VITREOUS BASE SHAVING.**

CUTTING RATE
5,000 cpm

**DOUBLE-
BLADE**
Long lasting,
excellent cutting
quality

**PERMANENTLY
OPEN**
Fast core vitrectomy

**CONSTANT
FLOW**
Controlled and
safe shaving

Geuder[®]
Precision made in Germany

UNO
COLORLINE
MACH2



THE NEW HIGH-SPEED VITREOUS CUTTER MACH2

IMPROVED SURGERY ROUTINE

Due to the permanently open cutting window, the GEUDER MACH2 vitreous cutter features a constant and high flow for a fast core vitrectomy.

MINIMIZED TRACTION

The flow, which is adjustable anytime, in combination with the maximum cutting rate of 5,000 cpm minimizes the tractions at the vitreous base. In connection with the reduced cutting window-to-tip distance, the MACH2 vitreous cutter allows a safe vitreous body shaving.

PREDICTABILITY = CONTROL

The stable, cutting rate-independent flow leads to a predictable behavior of the mobile retina, controllable by vacuum or flow.

PERFECT COMPLICATION MANAGEMENT

The constant and high flow provides an excellent holdability in cases of luxated lenses. The long-lasting cutting blades crush the nucleus safely even during long surgery times. Thanks to the constant and high flow, the minimized traction and improved control, the GEUDER MACH2 vitreous cutter is ideal in cases of traumata, bleedings into the vitreous cavity and retinal detachment.



UNO COLORLINE MACH2 23 G *

G-42302 - 1 pc per box

G-46302 - 6 pcs per box

*) fits Megatron/Megatron S3/Megatron S4

GEUDER AG reserves the right to make changes to technical details in response to recent developments. Geuder does not assume liability for the accuracy of each individual statement.

Illustrations not drawn to scale .