


OCU^{LED}

THE UNIQUE INTRAOCULAR
LED LIGHT SOURCE IN A HANDPIECE




DESIGNED FOR **megaTRON^{S4}**
PERFECTION THROUGH FLEXIBILITY.



- Sterile, single-use handpiece without an external illumination source
- Excellent lighting with optimal contrast
- Lower risk due to no IR and IV components
- Different models with spot, wide-angle and shielded designs available in 20, 23 and 25 Gauge

OCULED SINGLE-USE LED LIGHT SOURCES



Spot

Wide-angle

Wide-angle, shielded

20 Gauge

coolwhite

G-26401

G-26403

G-26405

neutralwhite

G-26402

G-26404

G-26406

23 Gauge

coolwhite

G-26431

G-26433

G-26435

neutralwhite

G-26432

G-26434

G-26436