

# DISPOSABLE ACCESSORY KITS

## FOR PHACOEMULSIFICATION



**ULTRASONIC TIP**  
30°, 1.17 MM

**G-24080  
DISPOSABLE PHACO ACCESSORY KIT**

Set consisting of:

- Ultrasonic tip; 30°, 1.17 mm, titanium
- Sleeve
- Phaco tip wrench and test chamber

10 pcs. per box, sterile

**BOTH SETS OFFER:**

- perfect balance of aspiration to irrigation
- easy handling
- hygienically safe, hence no risk of cross-contamination and infection
- economical due to no sterilization costs



**CMP  
ULTRASONIC TIP**  
30°, 0.80 MM

**G-24076  
DISPOSABLE CMP ACCESSORY KIT**

Set for Coaxial Micro Phaco consisting of:

- CMP tip, 30°, 0.8 mm, titanium
- CMP sleeve
- Phaco tip wrench and test chamber

10 pcs. per box, sterile

**ADDITIONAL ADVANTAGES OF CMP:**

- stable anterior chamber even at high vacuum up to 500 mmHg
- better surgical control and more freedom of movement
- less stress on ocular tissue
- rapid postoperative wound healing with less risk of leakage (due to 2.2 mm incision)
- reduction in surgically induced astigmatism